



Silica Reagent Set, 3-1000 µg/L SiO₂

Product #:

2581400

ZAR Price:

Contact Hach

Call for ship date

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

For ultralow range Silica determination by the Heteropoly Blue method.

Specifications

Instrument:	All except DR/3000, DR 900, DR/800 Series
Method:	PC II 8282
Number of tests:	40
Parameter:	Silica
Platform:	Solution
Range:	3 - 1000 µg/L SiO ₂
Sample size:	25 mL
Special Features:	UL, 40 tests

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)