



Ammonium p-toluenesulfonic acid, 25 g

| | |
|------------|--------------|
| Product #: | 2277924 |
| ZAR Price: | Contact Hach |
| Available | |

(Ammonium-PTSA) Protein-TKN standard: 7.402% N, 46.26% crude protein. For checking accuracy and digestion recovery in Total Nitrogen and TKN determinations. 25 g bottle. See 2277800 for the complete set of three standards with increasing resistance to digestion.

Specifications

Platform: Jar

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)