



## Nicotinic Acid p-Toluenesulfonic Acid, 25 g

Product #:	2278124
ZAR Price:	Contact Hach
Available	

(Nicotinic Acid PTSA) Protein-TKN standard: 4.743% N, 29.64% crude protein. For checking accuracy and digestion recovery in Total Nitrogen and TKN determinations. 25 g bottle. See 2277800 for the complete set of three standards with increasing resistance to digestion.

---

### Specifications

Platform: Jar

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)