



Primary standard, nitrogen, set of 3

Product #: 2277800
ZAR Price: Contact Hach
Ships within 2 weeks

Use this set for checking accuracy of crude protein tests and digestion recovery efficiency in Total Nitrogen and TKN methods. Set of three 25-g standards: Ammonium p-Toluenesulfonic Acid (Ammonium PTSA), Glycine PTSA and Nicotinic Acid PTSA.

Specifications

Description:	Primary Standard Set, for Kjeldahl Nitrogen Set of 3 Compounds: Ammonia p-Toluene Sulfonate, Glycine p-Toluene Sulfonate, Nicotinic Acid p-Toluene Sulfonate
Platform:	Jar
Volume:	25 g each

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)