



Glycine p-toluenesulfonic acid, 25 g

| | |
|---------------------|--------------|
| Product #: | 2278024 |
| ZAR Price: | Contact Hach |
| Ships within 3 days | |

(Glycine PTSA) Protein-TKN standard: 5.665% N, 35.40% crude protein. For checking accuracy and digestion recovery in Total Nitrogen and TKN determinations. 25 g bottle. See 2277800 for the complete set of three standards with increasing resistance to digestion.

Specifications

Platform: Jar

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)