



Cobalt Chloride Solution, 500 mL

Product #:

1422249

ZAR Price:

Contact Hach

Ships within 1 week

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Cobalt Chloride Solution, Bottle 500 mL.

For Zero point DO calibration for spectrophotometer linearisation.

Specifications

| | |
|---------------|--------------------------|
| Application: | Catalyst for Calibration |
| Description: | Catalyst |
| Formula: | CoCl_2 |
| Package Type: | Bottle |
| Platform: | Solution |
| Quantity: | 500 mL |