



Mixed Parameter Quality Control Standard (500 mL)

Product #: 2833149
ZAR Price: Contact Hach
Ships within 1 week

Quality Control Standard (waste water influent inorganics) for the parameters Ammonia, Nitrate, Phosphate, COD, Sulfate and TOC.

The standard solution contains the following concentrations:

- Ammonia: 15 mg/L $\text{NH}_3\text{-N}$
- Nitrate: 10 mg/L $\text{NO}_3\text{-N}$
- COD: 500 mg/L O_2
- Phosphate: 10 mg/L PO_4
- Sulphate: 400 mg/L SO_4
- TOC: 161 mg/L C

Specifications

Concentration 2:	Nitrate 10 mg/L $\text{NO}_3\text{-N}$
Concentration 3:	Phosphate 10 mg/L PO_4^{3-}
Concentration 4:	COD 500 mg/L O_2
Concentration 5:	Sulfate 400 mg/L SO_4^{2-}
Concentration 6:	TOC 161 mg/L C
Parameter:	Mixed
Platform:	Solution