



SpecCheck Gel secondary standard Kit, LR chlorine, DPD, 0-2.0 mg/L Cl_2

Product #:

2635300

ZAR Price:
Available

Contact Hach

Solid, ready for use Gel-Standards for instrument verification in the lab and in the field.

Ready for use SpecCheck Secondary Standards Kits include four 10 mL (1 inch) sealed glass cells with a blank and three different standards within the typical testing range, instructions, plastic case and Certificate of Analysis.

SpecCheck Standards are computer-formulated and have a two-year shelf life from date of manufacture.

Specifications

Concentration:	0 - 2.0 mg/L Cl_2
Instrument:	Includes values for all instruments
Method:	DPD
Parameter:	Chlorine
Quantity:	Set of 4 Vials