



## Ammonia Standard Solution, 10 mg/L, 500 mL

**Product #:** 15349  
**ZAR Price:** Contact Hach  
**Ships within 3 days**

### Ion Selective Electrode (ISE) standard solutions for laboratory analyses

ISE Standard Solution. Ammonia, 10.0 mg/L  $\pm 0.5$  mg/L  $\text{NH}_3\text{-N}$ . 500 mL bottle

#### Design for use with Hach meters portfolio and other meters.

The basis for accurate ISE measurements is good care and calibration of the Ion Selective Electrode.

#### Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

---

### Specifications

Concentration:	10.0 mg/L $\pm 0.5$ mg/L $\text{NH}_3\text{-N}$
Concentration 2:	0.00435 mol/L $\pm 0.00003$ mol/L $\text{Na}^+$
Package Type:	Bottle
Parameter:	Ammonium
Platform:	Solution
Quantity:	500 mL
Special Remarks:	Observed shelf life after opening, following recommended storage and use: 3 months
Storage conditions:	10 - 25 °C (protect from light)