



## Calcium Standard Solution, 10 mg/L, 100 mL

**Product #:** 2305442  
**ZAR Price:** Contact Hach  
**Ships within 3 days**

### Ion Selective Electrode (ISE) standard solutions for laboratory analyses

ISE Standard Solution. Calcium, 10.00 mg/L  $\pm 0.14$  mg/L  $\text{Ca}^{2+}$  (NIST). 100 mL bottle

#### Traceable to NIST standard reference materials

Always performing measurement results you can trust is essential. The quality of the standard solutions that are used are key.

#### Design for use with Hach meters portfolio and other meters.

The basis for accurate ISE measurements is good care and calibration of the Ion Selective Electrode.

#### Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

---

### Specifications

Concentration:	10.00 mg/L $\pm 0.14$ mg/L $\text{Ca}^{2+}$
Concentration 2:	25 mg/L $\pm 0.35$ mg/L $\text{CaCO}_3$
Package Type:	Bottle
Parameter:	Calcium
Platform:	Solution
Quantity:	100 mL
Special Remarks:	Observed shelf life after opening, following recommended storage and use: 3 months
Storage conditions:	10 - 25 °C (protect from light)