



Ammonia Standard Solution, 1000 mg/L, 1 L

Product #: 2354153
ZAR Price: Contact Hach
Ships within 3 days

Ion Selective Electrode (ISE) standard solutions for laboratory analyses

ISE Standard Solution. Ammonia, 1000 mg/L ± 20 mg/L $\text{NH}_3\text{-N}$. 1 L bottle

Design for use with Hach meters portfolio and other meters.

The basis for accurate ISE measurements is good care and calibration of the Ion Selective Electrode.

Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

Specifications

Concentration:	1000 mg/L ± 20 mg/L $\text{NH}_3\text{-N}$
Description:	Ammonia Standard Solution
Package Type:	Bottle
Parameter:	Ammonium
Platform:	Solution
Quantity:	1 L
Special Remarks:	Observed shelf life after opening, following recommended storage and use: 3 months
Storage conditions:	10 - 25 °C (protect from light)