



## Standard solution, sulfite, 15 mg/L as SO<sub>3</sub> (equivalent), 500 mL

**Product #:**

**2408449**

**ZAR Price:**

Contact Hach

**Ships within 3 days**

Sulfite equivalent standard solution. Prepared with stabilized Sodium Thiosulfate. Equivalent to 15 mg/L as SO<sub>3</sub>. 500 mL Bottle

---

### Specifications

Concentration:

15 mg/L

Description:

Sulfite Equivalent Solution as SO<sub>3</sub><sup>2-</sup> (NIST)

Platform:

Solution