



Quality control standard inorganics for drinking water, 500 mL

Product #: 2833049
ZAR Price: Contact Hach
Ships within 3 days

Inorganics Quality Control Standard for Drinking Water. Ready to use. No dilution required. Mixed standard includes: Fluoride: 1 mg/L F^- , Nitrate: 2 mg/L NO_3^- -N, Phosphate: 2 mg/L PO_4^{3-} , Sulfate: 50 mg/L SO_4^{2-} . 500 mL Bottle.

Specifications

Concentration:	Fluoride 1 mg/L F^-
Concentration 2:	Nitrate 2 mg/L NO_3^- -N
Concentration 3:	Phosphate 2 mg/L PO_4^{3-}
Concentration 4:	Sulfate 50 mg/L SO_4^{2-}
Description:	Drinking Water Inorganics, Mixed Parameter QC Standard
Platform:	Solution