



Quality control standard hardness, high range, 500 mL

Product #:

2833349

 **Hazardous**

ZAR Price:

Contact Hach

Ships within 3 days

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Ready to use - no dilutions required. Mixed standard includes: Total Hardness: 1000 mg/L as CaCO₃, Calcium Hardness: 500 mg/L as CaCO₃. 500 mL Bottle

Specifications

Concentration:

Total Hardness 1000 mg/L CaCO₃

Concentration 2:

Calcium Hardness 500 mg/L CaCO₃

Description:

High Range Hardness, Mixed Parameter QC Standard

Platform:

Solution