



Quality control standard hardness, low range, 500 mL

Product #: 2833449
ZAR Price: Contact Hach
Ships within 1 week

Ready to use - no dilutions required. Mixed standard includes: Total Hardness: 100 mg/L as CaCO_3 . Calcium Hardness: 50 mg/L as CaCO_3 .

Specifications

Concentration:	Total Hardness 100 mg/L CaCO_3
Concentration 2:	Calcium Hardness 50 mg/L CaCO_3
Description:	Low Range Hardness, Mixed Parameter QC Standard
Platform:	Solution