



## Quality control standard metals for drinking water, high range, 500 mL

**Product #:** 2833649  
**ZAR Price:** Contact Hach  
**Ships within 3-5 weeks**

Ready to use - no dilutions required. Mixed standard includes: Copper: 2.5 mg/L Cu, Iron: 1.5 mg/L Fe, Manganese: 5 mg/L Mn. 500 mL Bottle

---

### Specifications

Concentration:	Copper 2.5 mg/L Cu
Concentration 2:	Iron: 1.5 mg/L Fe
Concentration 3:	Manganese: 5.00 mg/L Mn
Description:	Metals Drinking Water (High Range) Mixed Parameter QC Standard
Platform:	Solution