



Chlorine Dosing Solution Ampules, 1000 - 1500 mg/L as Cl₂, 16/pk

Product #: 2504810
ZAR Price: Contact Hach
Ships within 1 week

Chlorine dosing solution for use with THM Formation Potential and Chlorine Demand Methods. 10 mL glass ampules. 16 ampules per package.

Specifications

Instrument:	DR 6000, DR 3900, DR 2800, DR 2700, DR 1900
Method:	10224
Method Number 2:	10223
Number of tests:	Varies
Platform:	Solution
Range:	1190 - 1310 mg/L as Cl ₂