



## Nitrate Standard Solution, 100 mg/L, 500 mL

**Product #:** 194749  
**ZAR Price:** Contact Hach  
**Available**

### Ion Selective Electrode (ISE) standard solutions for laboratory analyses

ISE Standard Solution. Nitrate, 100 mg/L  $\pm 1$  mg/L  $\text{NO}_3\text{-N}$  (NIST). 500 mL bottle

#### Traceable to NIST standard reference materials

Always performing measurement results you can trust is essential. The quality of the standard solutions that are used are key.

#### Design for use with Hach meters portfolio and other meters.

The basis for accurate ISE measurements is good care and calibration of the Ion Selective Electrode.

#### Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

---

### Specifications

Concentration:	100 mg/L $\pm 1$ mg/L $\text{NO}_3\text{-N}$
Concentration 2:	443 $\pm 4$ mg/L $\text{NO}_3$
Description:	Nitrogen-Nitrate Standard Solution as $\text{NO}_3\text{-N}$ (NIST)
Package Type:	Bottle
Parameter:	Nitrate
Platform:	Solution
Quantity:	500 mL
Special Remarks:	Observed shelf life after opening, following recommended storage and use: 3 months
Storage conditions:	10 - 25 °C (protect from light)