



Calcium Standard Solution, 100 mg/L, 500 mL

Product #: 4457649
ZAR Price: Contact Hach
Ships within 1 week

Ion Selective Electrode (ISE) standard solutions for laboratory analyses

ISE Standard Solution. Calcium, 100 mg/L ± 1 mg/L Ca^{2+} (NIST). 500 mL bottle

Traceable to NIST standard reference materials

Always performing measurement results you can trust is essential. The quality of the standard solutions that are used are key.

Design for use with Hach meters portfolio and other meters.

The basis for accurate ISE measurements is good care and calibration of the Ion Selective Electrode.

Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

Specifications

Concentration:	100 mg/L ± 1 mg/L Ca^{2+}
Package Type:	Bottle
Parameter:	Calcium
Platform:	Solution
Quantity:	500 mL
Special Remarks:	Observed shelf life after opening, following recommended storage and use: 3 months
Storage conditions:	10 - 25 °C (protect from light)