



## AluVer 3 Aluminum Reagent Powder Pillows

<b>Product #:</b>	<b>1429099</b>
ZAR Price:	Contact Hach
Ships within 3 days	

Powder reagent in sealed package for the determination of Aluminium, 0.008-0.800mg/L Al, Aluminon Method. Pack of 100.

---

### Specifications

Method:	Aluminon
Number of tests:	100
Platform:	Powder Pillows
Quantity:	100 /pk

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)