



Aluminium Reagent Set, Powder Pillows, 0.008-0.800mg/L Al

Product #:

2242000

ZAR Price:
Available

Contact Hach

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Powder reagents in sealed packages, for determination of aluminium. 100 tests.

Specifications

Description:	Aluminum Reagent Set, Powder Pillows, Aluminon
Instrument:	All
Method:	8012
Number of tests:	100
Parameter:	Aluminum
Platform:	Powder Pillows
Range:	0.008 - 0.800 mg/L Al
Special Remarks:	Set contains: AluVer 3 Aluminium Powder Pillows (1429099), Ascorbic Acid Powder Pillows (1457799), Bleaching 3 Powder Pillows (1429449)
Storage conditions:	10 °C to 30 °C

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)