



Aluminum Reagent Set, 0.002-0.250 mg/L Al

Product #:

2603700

ZAR Price:

Contact Hach

Ships within 3 days

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

For the determination of Aluminium, Eriochrome Cyanine R (ECR) method. For 10 mL samples. 100 tests.

Specifications

Instrument:	DR6000, DR3900, DR1900
Method:	8326
Number of tests:	100
Parameter:	Aluminum
Range:	0.002 - 0.250 mg/L Al

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)