



ECR Reagent Powder Pillows, pk/100

Product #: 2603849
ZAR Price: Contact Hach
Ships within 3 days

Powder reagent in sealed package for the determination of Aluminium, Eriochrome Cyanine R Method. For 20mL sample.

Specifications

Method:	Eriochrome Cyanine R
Number of tests:	100
Parameter:	Aluminum
Platform:	Powder Pillows
Quantity:	100/pk

Required Accessories

- Hexamethylenetetramine buffer reagent, 100 powder pillows (Item 2603999)
- ECR Masking Reagent Solution (25mL SCDB) (Item 2380123)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)