



## ECR Masking Reagent Solution (25mL SCDB)

**Product #:** 2380123  
**ZAR Price:** Contact Hach  
**Ships within 3 days**

For the determination of Aluminium, Eriochrome Cyanine R Method. Self-contained dropping bottle.

---

### Specifications

Method:	Eriochrome Cyanine R
Number of tests:	100
Platform:	Solution
Range:	0.002 - 0.250 Al

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)