



## Solution, lead acetate, 10%, 100 mL

**Product #:** 1458042  
**ZAR Price:** Contact Hach  
**Ships within 1 week**

Reagent for arsenic and selenium determination. APHA formulation

---

### Specifications

Instrument:	All except DR/800, PC II
Method:	8013
Number of tests:	Varies
Platform:	Solution
Quantity:	100
Range:	0 - 0.200 mg/L As

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)