



Barium Reagent Powder Pillows, 2-100mg/L Ba

Product #: 1206499
ZAR Price: Contact Hach
Ships within 3 days

Powder reagent in sealed packages for the determination of Barium, turbidimetric method. 100/pk

Specifications

Instrument:	All except PC II
Method:	Turbidimetric
Number of tests:	100
Parameter:	Barium
Platform:	Powder Pillows
Quantity:	100/pk
Range:	2 - 100 mg/L Ba

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)