



## Reagent powder pillows, ChromaVer 3, chromium , 25 mL, pk/1000

Product #: 1206628

 **Hazardous**

**OBSOLETE ITEM**  
This item is no longer available.

Items with this mark may be considered hazardous under some shipping conditions.

Suggested replacements:  
• 1206699

If necessary, we will change your selected shipping method to accomodate these items.

For hexavalent chromium determination by the Diphenylcarbohydrazide Method. USEPA accepted for reporting for wastewater analyses using Hach Method 8023. Range: up to 0.70 mg/L Cr. Sample Size: 25 mL.

---

### Specifications

EPA compliant:	Yes
Method:	1,5-Diphenylcarbohydrazide
Platform:	Powder Pillows

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)