



Reagent powder pillows, ChromaVer 3, chromium, 5 or 10 mL, pk/1000

Product #:

1271028

ZAR Price:
Available

Contact Hach

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

For determination of hexavalent chromium by the Diphenylcarbohydrazide method. USEPA accepted for reporting for wastewater analyses using Hach Method 8023, Diphenylcarbohydrazide Method. Photometric Sample Size: 10 mL. Photometric Range: up to 0.70 mg/L. Kit Sample Size: 5 mL. Kit Range: Up to 1.50 Pack of 1000 powder pillows.

Specifications

EPA compliant:	Yes
Instrument:	All except DR/3000, DR/2000, DR 900
Method:	8023
Number of tests:	1000
Parameter:	Chromium
Platform:	Powder Pillows
Range:	0.010 - 0.700 mg/L Cr ⁶⁺

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)