



Digital titrator cartridge, sodium thiosulfate, 0.2068N

Product #: 2267601
ZAR Price: Contact Hach
Available

For determination of Chromate by Digital Titrator titration. 0.2068 N $\text{Na}_2\text{S}_2\text{O}_3$. Digits x 0.2 = mg/L CrO_4 in a 50 mL sample.

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specifications

Method:	Iodometric
Number of tests:	40 - 100
Parameter:	Chromate
Platform:	Digital Titrator
Range:	20 - 400 mg/L CrO_4
Titant:	0.2068 N $\text{Na}_2\text{S}_2\text{O}_3$

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)