



Reagent set, chromate, digital titrator

Product #: 2272400

 **Hazardous**

OBSOLETE ITEM
This item is no longer available.

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

For determination of chromate by sodium thiosulfate Digital Titrator titration. Hach Method 8211. Range: 20-400 mg/L CrO₄. Reagent Set includes: Dissolved Oxygen 3 Reagent Powder Pillows, Potassium Iodide Powder Pillows, Sodium Thiosulfate Cartridge (0.2068 N) and Starch Indicator Solution. Approximately 40-100 tests

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specifications

Method:	8211
Number of tests:	40 - 100
Range:	20 - 400 mg/L CrO ₄ ²⁻

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)