



Be Right™



Standard solution, sodium thiosulfate, stabilized, 0.0113 N, 100 mL MDB

Product #:

2409032

ZAR Price:
Available

Contact Hach

Stabilized sodium thiosulfate titrant for chromate determination

Specifications

Platform:

Solution

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)