



Porphyrim 2 Reagent Powder Pillows, 25 mL, pk/100

Product #:

2187569

ZAR Price:

Contact Hach

Ships within 3-5 weeks

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Powder reagent in sealed package for the determination of Copper by the Porphyrim method, 2 - 210 µg/L Cu. For 25 mL sample.

Specifications

Method:	Porphyrim
Platform:	Powder Pillows
Quantity:	100 /pk

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)