



## Reagent set, copper, free & total

**Product #:** 2439200  
**ZAR Price:** Contact Hach  
**Ships within 3 days**

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

For determination of free and complexed (chelated) Copper by the Binchoninate method. Range: up to 5.00 mg/L. Set includes: Free Copper Reagent Powder Pillows and Hydrosulfite Powder Pillows for 100 tests.

---

### Specifications

Instrument:	DR 6000, DR 3900, DR 1900, DR 900, PCII
Method:	Bicinchoninate
Number of tests:	100
Platform:	Powder Pillows
Range:	0.04 - 5.00 mg/L Cu

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)