



Ferrous Iron Reagent Powder Pillows, 0.02-3.00mg/L Fe (II)

Product #: 103769
ZAR Price: Contact Hach
Ships within 3 days

Powder reagent in sealed packages for the determination of Iron in the ferrous form (Fe^{2+}). pk/100.

Specifications

Instrument:	All except PC II
Method:	1, 10 Phenanthroline
Number of tests:	100
Parameter:	Iron, Ferrous
Platform:	Powder Pillows
Quantity:	100 /pk
Range:	0.02 - 3.00 mg/L Fe^{2+}

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)