



## FerroVer Iron Reagent Powder Pillows (pk/100)

<b>Product #:</b>	<b>92799</b>
ZAR Price:	Contact Hach
Available	

Powder reagent in sealed packages for the determination of Total Iron, 1,10 Phenanthroline Method. For 5 mL sample.

---

### Specifications

Method:	FerroVer
Parameter:	Iron
Platform:	Powder Pillows
Quantity:	100

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)