



FerroZine Iron Reagent Solution Pillows, 0.009-1.400mg/L Fe

Product #:

230166

ZAR Price:

Contact Hach

Ships within 3 days



Hazardous

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

For the determination of soluble ferric (Fe^{3+}) and ferrous (Fe^{2+}) Iron at trace levels.

FerroZine Iron Reagent determines soluble ferric (Fe^{3+}) and ferrous (Fe^{2+}) Iron. Digestion is needed for strongly chelated Iron complexes (such as those with EDTA) and Total Iron. For trace level Iron analysis. 50/pack

Specifications

Instrument:	All except DR 900, PC II
Method:	8147
Number of tests:	50
Parameter:	Iron
Quantity:	50 /pk
Range:	0.009 - 1.400 mg/L Fe

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)