



FerroZine Iron Reagent Solution, 0.009-1.400mg/L Fe

Product #:

230149

ZAR Price:
Available

Contact Hach

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

For the determination of soluble ferric (Fe^{3+}) and ferrous (Fe^{2+}) Iron at trace levels.

Specifications

Instrument:	All except DR 900, PC II
Method:	8147
Number of tests:	500 - 1000
Parameter:	Iron
Platform:	Solution
Quantity:	500 mL
Range:	0.009 - 1.400 mg/L Fe

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)