



## Dissolved Iron Chemkey Reagents (box of 25)

<b>Product #:</b>	<b>8636000</b>
ZAR Price:	Contact Hach
Ships within 1 week	

Dissolved Iron Chemkeys, range 0.05 - 3.00 mg/L Fe. For SL1000 Portable Parallel Analyser (PPA).

---

### Specifications

Parameter:	Dissolved Iron
Platform:	Chemkeys
Quantity:	pk/25
Range:	0.05 - 3.00 mg/L Fe
What's included?:	25 Chemkeys

---

### What's included?

25 Chemkeys