



FerroVer Iron Reagent, Accuvac, 0.02-3.00 mg/L Fe

Product #:	2507025
ZAR Price:	Contact Hach
Ships within 3 days	

Accuvac ampules contain the precise amount of reagent for a single test and can be used as a measurement cuvette. For determination of total dissolved Iron, as well as complexed iron compounds, by the 1,10 Phenanthroline Method.

Specifications

EPA compliant:	Yes
Method:	8008
Number of tests:	25
Parameter:	Iron
Platform:	Accuvac
Range:	0.02 - 3.00 mg/L Fe

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)