



## Lead reagent set, 5-150 µg/L Pb

**Product #:**

**2375000**

**ZAR Price:**

Contact Hach

**Ships within 3 days**

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Acid soluble Lead in a potable water sample is first concentrated on a Fast Column Extractor. The Lead is then eluted from the Extractor and determined colorimetrically with an indicator.

---

### Specifications

Instrument:	All except DR 900
	DR/800 Series
Method:	8317
Number of tests:	20
Parameter:	Lead
Range:	5 - 150 µg/L Pb
Shipping Method:	Lead Reagent set

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)