



## Manganese Reagent Set, Powder Pillows, periondate oxidation, 0.1-20.0 mg/L Mn, 10 mL sample

Product #:

2430000

ZAR Price:

Contact Hach

Ships within 3 days

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Powder reagents in sealed packages for the determination of soluble Manganese, Periodate Oxidation Method. For 10mL sample. 100 tests

---

### Specifications

EPA compliant:	Yes
Instrument:	All except DR/3000, DR/2000
Method:	8034
Number of tests:	100
Parameter:	Manganese
Platform:	Powder Pillows
Range:	0.1 - 20.0 mg/L Mn
Sample Volume:	10 mL