



Manganese Reagent Set, 0.006-0.700 mg/L Mn, 25 mL

Product #:

2243300

ZAR Price:

Contact Hach

Ships within 3 days

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

For the determination of Manganese with PAN Method. For 25 mL sample.

For determination of Manganese by the PAN Method. Hach Method 8149. Range: 0.006 - 0.700 mg/L Mn. Sample Size: 25 mL. Reagent Set includes: Ascorbic Acid Powder Pillows, Alkaline-Cyanide Reagent and PAN Indicator Solution (0.1%) for 100 tests.

(DR/3000, DR/2000)

Specifications

Instrument:	DR/3000, DR/2000
Method:	8149
Number of tests:	100
Range:	0.006 - 0.700 mg/L Mn
Sample Volume:	25 mL
Special Features:	Manganese Reagent set, 25ml