



Manganese Reagent Set, PAN, 0.006-0.700mg/L Mn

Product #:

2651700

ZAR Price:

Contact Hach

Ships within 3 days

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Self contained dropping bottles and powder reagents in sealed packages for determination of Manganese by the PAN method.

Set contains: Alkal. Cyanide reagent (2122326) 50mL, Ascorbic acid (1457799), PAN-Indicator solution (2122426) 50 mL

Specifications

Instrument:	All except DR/3000, DR/2000, DR/850, DR/820
Method:	PAN
Number of tests:	50
Parameter:	Manganese
Range:	0.006 - 0.700 mg/L Mn
Sample Volume:	50 mL
Special Remarks:	Set contains: Alkal. Cyanide reagent (2122326) 50 ml, Ascorbic acid (1457799), PAN-Indicator solution (2122426) 50 ml