



Stannous sulfate solution, HgEx reagent A, 20 mL ampules, pk/25

Product #:

2658825

 **Hazardous**

OBSOLETE ITEM
This item is no longer available.

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

For mercury determination by the Cold Vapor Mercury Concentration Method. Hach Method 10065. Reagent is a component of reagent set 2658300.

Specifications

Method:	Cold Vapor Concentration
Number of tests:	25
Platform:	Solution

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)