



Be Right™



Sodium hypochlorite solution, HgEx reagent C, 55 mL

Product #: 2659059
ZAR Price: Contact Hach
Available

For mercury determination by the Cold Vapor Mercury Concentration Method. Hach Method 10065. Reagent is a component of reagent set 2658300.

Specifications

Method: Cold Vapor Concentration
Number of tests: 25
Platform: Solution

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)