



Boron Reagent Powder Pillows, 0.2-14.0mg/L B

Product #: 1417099
ZAR Price: Contact Hach
Ships within 3 days

Powder reagent in sealed packages for the determination of Boron with Carmine Method. pk/100.

Specifications

Instrument:	All except DR/800, PC II
Method:	Carmine
Number of tests:	100
Parameter:	Boron
Quantity:	100
Range:	0.2 - 14.0 mg/L B

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)