



**Be Right™**



## **Boron Reagent Powder Pillows, 0.2-14.0mg/L B**

**Product #:** 1417099  
**ZAR Price:** Contact Hach  
**Available**

Powder reagent in sealed packages for the determination of Boron with Carmine Method. pk/100.

---

### **Specifications**

Instrument:	All except DR/800, PC II
Method:	Carmine
Number of tests:	100
Parameter:	Boron
Quantity:	100
Range:	0.2 - 14.0 mg/L B

---

### **Required Accessories**

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)