



Be Right™



MolyVer 1 Molybdenum Powder Pillows

Product #: 1414669
ZAR Price: Contact Hach
Available

Powder reagent in sealed package for the determination of high range molybdenum (0-40.0mg/L Mo), Mercaptoacetic Acid Method. For 25mL sample.

Specifications

| | |
|------------------|----------------|
| Method: | 8036 |
| Number of tests: | 100 |
| Parameter: | Molybdenum |
| Platform: | Powder Pillows |
| Quantity: | 100/pk |

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)