



MolyVer 2 Molybdenum Powder Pillows

Product #:	1414869
ZAR Price:	Contact Hach
Ships within 1 week	

Powder reagent in sealed package for the determination of high range molybdenum (0-40.0mg/L Mo), Mercaptoacetic Acid Method. For 25mL sample. Pack of 100.

Specifications

Method:	8036
Number of tests:	100
Parameter:	Molybdenum
Platform:	Powder Pillows
Quantity:	100/pk

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)