



Molybdenum Reagent Set, Powder Pillows, 0-35.0 mg/L Mo

Product #: 2243400
ZAR Price: Contact Hach
Ships within 3 days

Powder reagents in sealed packages for the determination of Molybdenum, Colorimetric Mercaptoacetic Acid Method. For 25 mL sample.

For determination of Molybdenum by the Mercaptoacetic Acid Method (Hach method 8036). Range: Up to 35.0 mg/L as Mo_6^+ . Sample Size: 25 mL. Set Contains: MolyVer 1, 2 and 3 Molybdenum Reagent Powder Pillows for 100 tests.

Specifications

Instrument:	DR/3000, DR/2000
Method:	8036
Number of tests:	100
Platform:	Powder Pillows
Range:	0 - 35.0 mg/L Mo
Sample size:	25 mL

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)